


CyberMed Rx

Desktop Docking Station



Technical Specifications

I/O Ports	2x USB 2.0 ports 1x RS232 Serial port 1x HDMI – Out 1x Gigabit LAN port 1x Power Jack
Theft Deterrent	Kensington Lock Support
Relative Humidity	10%~90% (non-condensing)
Dimensions	9.7" x 10.4" x 8.3" (L, W, D)
Weight	6 lbs without the tablet or additional batteries, 6.5lbs with two batteries
Environment	Ambient Temperature: 0°C~40°C (operating)
Certifications	FCC Class B, CE, RoHS
Battery Capacity	2 slots + 2 in tablet for a total of 4 batteries simultaneously charging
Battery Recharge Time	2 hours to recharge all batteries simultaneously